# Return-to-Walk A

## Fußheberorthese

AFO

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## User Manual

[www.wjt-ortho.com](http://www.wjt-ortho.com)

ORTHO

## ◆ Inhalt

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Achtung

## ◆ Content

Description	EN-2
Indication	EN-2
Contraindication	EN-2
Function	EN-2
Reimbursement Code	EN-2
Medical Device	EN-2
Safety	EN-3
Specific Safety Instructions	EN-4
REF/Sizes	EN-4
Scope of Delivery	EN-4
Application/Adjustment	EN-5
Comfort/Hygiene	EN-6
Material	EN-6
Produkt Care/Cleaning/Maintenance	EN-6
Storage/Disposal	EN-6
CE-Conformity	EN-6
Legend to Symbols	EN-6

## ◆ Description

Return-to-Walk A is an AFO made from Carbon. It supports patients with weak dorsiflexion due to a neurological or neuro-orthopedic disease or due to a traumatic event. Especially patients with an additional weakness of the quadriceps will profit from the anterior tibia frame of the orthosis. The patient should exhibit only light spasticity. The orthosis should be used as an orthopedic aid for the lower extremities only. It should not be worn directly on the skin but rather over a thin stocking/sock.

## ◆ Indication

- After a stroke
- Neural muscular atrophy
- Poliomyelitis
- Peroneal lesion
- Multiple sclerosis
- Traumatic brain injury

## ◆ Contraindication

- Distinct foot supination, caused by e.g. severe spasticity
- Foot deformities
- Moderate to severe edema formation in the lower extremities

## ◆ Function

Prevention of foot drop during swing phase with simultaneous support of torque during knee extension, to minimize the risk of falling. In addition, the patient will regain a physiological gait pattern and a quicker pace.

## ◆ Reimbursement Code

DE HMV-Nr.: 23.03.02.6027

AT Tarif-Nr.: 34382 038

## ◆ Medical Device

The Return-to-Walk A AFO is a medical device.

## ◆ ⚠ Safety

- Before this device is used, the orthosis has to be individually adapted to the patient by an experienced orthopedic technician and the patient has to be trained in the orthosis' use according to the doctor's instructions.
- Incorrect selection, adaptation, application, and/or failure to monitor the correct fit of this product may cause health damage.
- Check orthosis for damage (such as loose connections, fasteners, belts, etc.) each time it is put on.
- Do not remove the textile label from textile products.
- In case you suspect any performance malfunctions (e.g. an apparent reduction or lack of effect), immediately call on your doctor or orthopedic technician.
- All serious incidents that occur in connection with this product have to be reported to the Wilhelm Julius Teufel GmbH and the competent authority of the member state in which you are established. A „serious incident“ is any incident that caused, may have caused, or may cause, directly or indirectly, one of the following consequences:
  - a) death of a patient, user, or other persons
  - b) temporary or permanent health deterioration of a patient, user or other persons
  - c) a serious public health threat.
- Never wear the orthosis on broken skin.
- The orthosis must be checked for proper function, correct fit and potential damage by a doctor or orthopedic technician at regular intervals. The timing of these intervals has to be defined individually for each patient.
- Should the product have been subject to disproportionate stress, it must be checked for potential damage by an orthopedic technician before further use.
- Harmful pressure caused by the orthosis may cause skin injuries or disrupt blood flow and must be avoided.
-  This product is intended for use on one person. It may not be reused on another person.

## ◆ ⚠ Specific Safety Instructions

- Avoid local pressure points and possible concurrent disruption of blood flow by not fastening the orthosis too tightly.
- The orthosis can only be used for a limited range of athletic activities. If the patient wants to perform any athletic activities, the orthopedic technician should inform her or him about the loading capacity and resilience of the orthosis.
- Movements during which the orthosis gets bent excessively (e.g. initial contact of foot with the whole forefoot instead of the tip of the toes while climbing/descending stairs, or squatting) should be avoided.
- Maximum body weight is 100 kg.

## ◆ REF/Sizes

Return-to-Walk A				
REF Left	REF Right	Size	Shoe Size	Height of Orthosis
22.390.312.00	22.390.412.00	S	35-38	35 cm
22.390.313.00	22.390.413.00	M	38-41	36 cm
22.390.314.00	22.390.414.00	L	41-44	37 cm
22.390.315.00	22.390.415.00	XL	44-47	37 cm
Optional Accessories				
REF	Product			
22.390.950.00	Return-to-Walk A Pad&Strap Kit S-M			
22.390.950.01	Return-to-Walk A Pad&Strap Kit L-XL			

## ◆ Scope of Delivery

Carbon orthosis with tibia pad. Two straps, which can be closed with a hook-and-loop closure, attach the orthosis to the lower leg. An insole with double-sided adhesive tape to attach the sole and a user manual with all relevant information are included, to ensure the ideal adjustment of the orthosis to the patient.

## ◆ Application/Adjustment

1. Determine the shoe size of the patient and select the size of the orthosis accordingly.
2. For best results the patient should be wearing a sturdy outdoor shoe with a strong heel counter. Heel height should be about 1 cm.
3. Take the sole out of the outdoor shoe, and cut the foot plate of the orthosis to the right size (Fig. 1).

4. To achieve a very good fit, the width of the orthosis' foot plate should only be reduced very carefully. The length of the foot plate depends on the position of the lateral strut, which should be positioned at the height of the metatarsal-5 base. There should be not contact with the ankle joint. The core of the sole is made of carbon, hence only the edges of the sole should be trimmed. (Fig. 2).

5. To attach the orthosis to the patient's lower leg, open both straps. Place the foot of the patient in the right position, so the tibia touches the tibia pad. Close the upper strap, then the lower strap. The straps' length can be adjusted by cutting them, and changing the position of the hook-and-loop closure accordingly (Fig. 3).

## ◆ Comfort/Hygiene

Do not wear the orthosis directly on the skin, but instead over thin clothing (e.g. a sock of appropriate length).

## ◆ Material

The orthosis is made from carbon fiber composite material. The textile pad is made of soft fleece with air pockets. The straps are made of soft velours with hook-and-loop closures. The product is latex-free. The insole is basically foam material.

## ◆ Produkt Care/ Cleaning/Maintenance

### • Textile Component:

Hand wash only.

Do not bleach.

Do not tumble dry.

Do not iron.

Do not dry-clean.

### • Carbon Fiber Composite Component:

Clean with a wet cloth.

## ◆ Storage/Disposal

Store the product in a dry place. Protect from dust, dirt and moisture. This product does not contain any environmentally hazardous substances. It can be disposed of with household waste.

## ◆ CE-Conformity

This medical device complies with the requirements of the MDR, Regulation (EU) 2017/745.

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