Ankle Air PRO

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User Manual

Inhalt

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Contents

Description	EN-2
Indication	EN-2
Contraindication	EN-2
Function	EN-3
Reimbursement Code	EN-3
Medical Device	EN-3
Safety	EN-3
REF/Sizes	EN-4
Scope of Delivery	EN-4
Application/Adjustment	EN-4
Comfort/Hygiene	EN-5
Materials	EN-5
Product Care/Cleaning/Maintenance	EN-6
Storage/Disposal	EN-6
Legend to Symbols	EN-6

Description

The Ankle Air PRO ankle brace is an orthosis which may be used both during the post-operative stage and for injury prevention, and helps users prevent twisting their foot. Ligaments, soft tissue and osseus structures of the lower and upper ankle joint are stabilized, which furthers the healing process of the injured structures.

The orthosis can be used during everyday life or when doing sports to prevent positions and movements during which people typically twist their ankles. This way, the orthosis helps prevent injuries of the ankle joint.

Because the orthosis can be thermoformed, the novel, patented shape of the brace snugly fits the patients' individual ankle anatomies.

The innovative straps are very thin, and can be positioned both easily and quickly. The brace does not add any bulk, and is hardly noticeable when worn.

The degree of joint stabilization can be adjusted individually. By changing the adjustment of the brace, the mobility of the ankle can be increased, which makes the brace downgradable. The textile on the central stabilizing element is highly breathable, and prevents the accumulation of heat and moisture.

◆ Indication

- Postoperative or posttraumatic treatment of the ligamentous apparatus of the ankle joint
- Acute injuries to the ligaments of the upper and lower ankle joint
- Chronic instability of the upper and lower ankle joint
- Distorsions

Contraindication

- General lesions of the skin, redness of the skin where the brace will be worn
- Severely disturbed blood flow of the lower extremities
- Lymphatic problems on the affected foot or leg

Function

Securing the ankle joint against unphysiological movements which may cause an injury through twisting the lower or upper ankle. This function is fulfilled both during the posttraumatic/postoperative period and for injury prevention.

Reimbursement Code

DE: HMV: 23.02.02.2013

Medical Device

The Ankle Air PRO is a medical device.

◆ **A** Safety

- Before this device is used, the orthosis
 has to be individually adapted to the
 patient by an experienced orthopedic
 technician and the patient has to be
 trained in the orthosis' use according to
 the doctor's instructions.
- Incorrect selection, adaptation, application, and/or failure to monitor the correct fit of this product may cause health damage.
- Check orthosis for damage (such as loose connections, fasteners, belts, etc.) each time it is put on.
- Do not remove the textile label from textile products.
- In case you suspect any performance malfunctions (e.g. an apparent reduction or lack of effect), immediately call on your doctor or orthopedic technician.

- All serious incidents that occur in connection with this product have to be reported to the Wilhelm Julius Teufel GmbH and the competent authority of the member state in which you are established. A "serious incident" is any incident that caused, may have caused, or may cause, directly or indirectly, one of the following consequences:
 - a) the death of a patient, user, or other persons
 - b) temporary or permanent health deterioration of a patient, user or other persons
 - c) a serious public health threat.
- Never wear the orthosis on broken skin.
- The orthosis must be checked for proper function, correct fit and potential damage by a doctor or orthopedic technician at regular intervals. The timing of these intervals has to be defined individually for each patient.
- Should the product have been subject to disproportionate stress, it must be checked for potential damage by an orthopedic technician before further use.
- Harmful pressure caused by the orthosis may cause skin injuries or disrupt blood flow and must be avoided.
- This product is intended for use on one person. It may not be reused on another person.

REF/Sizes

REF Left	REF Right	Size	Body height
22.470.000.02	22.470.000.12	S-M	< 165 cm
22.470.000.03	22.470.000.13	L-XL	> 165 cm


Scope of Delivery

Ankle brace consisting of three straps with hook-and-loop closures and a user manual

◆ Application/Adjustment

- **1.** First, determine the patient's height to select the right brace size for the patient.
- 2. Place the brace against the medial malleolus. The position of the ankle displays an angle of 90°.

To acchieve an optimum fit of the brace, the central stabilizing element, which consists of fiber composite material, can be thermoformed (fig. 1a).

4. Place strap no. 2 around the back of the ankle, and attach it to the front on the medial side with the two hook dots that have already been installed there (fig. 2). The strap can be adjusted in length by cutting it (fig. 1b).

5. Wrap strap no. 3 once around the bottom of the foot (fig. 3a).

Now place the strap across the back of the foot, then one more time around the bottom of the foot. The anatomically shaped loop strap element (A) should be placed directly onto the lateral side of the foot to prevent the brace from shifting out of place (fig. 3b).

Finally, place the strap diagonally across the hook-and-loop strap no. 2, and attach it on the hook-and-loop strap no.1 (fig. 3c).

Where the straps cross, loop strap no. 3 should be attached to the hook dot on strap no. 2, to continuously ensure that the strap pulls in the correct direction. The length of the strap both in front and behind the anatomically shaped lopp strap element (A) can be adjusted by cutting the strap to the desired length. The orthosis can be downgraded and the range of motion, especially of the lower ankle joint, can be increased by removing the loop strap element (A) and the strap following it, and changing the way the straps wrap around the ankle and foot. This adjustment should be administered by an orthopedic technician. The correct time to downgrade the orthosis should be determined by the attending physician.

Comfort/Hygiene

Do not wear the orthosis directly on the skin, but instead over thin clothing (e.g. a sock of appropriate length).

Materials

Fiber composite
Polyethylene
Spacer fabric
Velours
Hook-and-loop material

The orthosis is latex-free.

Product Care/Cleaning/ Maintenance

Hand wash only.

Do not bleach.

Do not tumble dry.

Do not iron.

Do not dry clean.

Storage/Disposal

Store the product in its original packaging in a dry place. Protect from dust, dirt and moisture. This product does not contain any environmentally hazardous substances. It can be disposed of with regular household waste.

◆ Legend to Symbols

Manufacturer

Date of Manufacture

Distributor

Importer

EC Authorized Representative

Lot Number

Part Number

Serial Number

UDI Data Matrix

Follow User Manual

Caution

CE-Marking According to European Union Council Regulation 2017/745 Concerning Medical Devices (MDR)

Single Patient, for Multiple Use

Store in a Dry Place, Protect from Moisture

Protect from Sunlight

Temperature Limit

Humidity Limit

Atmospheric Pressure Limit, for Transport and Storage

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Wilhelm Julius Teufel GmbH Robert-Bosch-Straße 15 73117 Wangen Deutschland/Germany Phone: +49 (0)7161 15684-0