
Return-to-Walk P & P/Plus Fußheberorthese

Gebrauchsanweisung User Manual

◆ Inhalt

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Content

Description	EN-2
Indication	EN-2
Contraindication	EN-2
Function	EN-2
Reimbursement Code	EN-2
Medical Device	EN-2
Safety	EN-3
Specific Safety Instructions	EN-4
REF/Sizes	EN-4
Scope of Delivery	EN-5
Application/Adjustment	EN-5
Comfort/Hygiene	EN-6
Material	EN-6
Produkt Care/Cleaning/Maintenance	EN-6
Storage/Disposal	EN-6
CE-Conformity	EN-6

Description

Return-to-Walk P and P/Plus are AFOs made from carbon. They support patients with weak dorsiflexion due to a neurological or neuro-orthopedic disease or due to a traumatic event. The patient should be only slightly spastic. The orthoses should be used as an orthopedic aid for the lower extremities only. They should not be worn directly on the skin, but always over a thin sock or stocking.

The Return-to-Walk P/Plus is more rigid than the Return-to-Walk P. The patient and the orthopedic technician should choose the adequate version for the patient together.

Indication

- Moderate drop foot
- Peroneal paralysis without spasticity
- Peroneal paralysis with spasticity
- Follow-up treatment of stroke
- Hemiplegia
- Poliomyelitis

Contraindication

- Distinct foot supination, caused e.g. by severe spasticity
- Foot deformities
- Moderate to severe edema formation in the lower extremities

Function

Prevention of foot drop during swing phase to minimize the risk of falling. In addition, the patient will regain a physiological gait pattern and a quicker pace.

Reimbursement Code

DE HMV-Nr.: 23.03.02.6025 **AT** Tarif-Nr.: 34382 039

Medical Device

The Return-to-Walk P (/Plus) is a medical device.

Safety

- Before this device is used, the orthosis has to be individually adapted to the patient by an experienced orthopedic technician and the patient has to be trained in the orthosis' use according to the doctor's instructions.
- Incorrect selection, adaptation, application, and/or failure to monitor the correct fit of this product may cause health damage.
- Check orthosis for damage (such as loose connections, fasteners, belts, etc.) each time it is put on.
- Do not remove the textile label from textile products.
- In case you suspect any performance malfunctions (e.g. an apparent reduction or lack of effect), immediately call on your doctor or orthopedic technician.
- All serious incidents that occur in connection with this product have to be reported to the Wilhelm Julius Teufel GmbH and the competent authority of the member state in which you are established. A "serious incident" is any incident that caused, may have caused, or may cause, directly or indirectly, one of the following consequences:
 - a) the death of a patient, user, or other persons
 - b) temporary or permanent health deterioration of a patient, user or other persons
 - c) a serious public health threat.
- Never wear the orthosis on broken skin.
- The orthosis must be checked for proper function, correct fit and potential damage by a doctor or orthopedic technician at regular intervals. The timing of these intervals has to be defined individually

- for each patient.
- Should the product have been subject to disproportionate stress. it must be checked for potential damage by an orthopedic technician before further use.
- Harmful pressure caused by the orthosis may cause skin injuries or disrupt blood flow and must be avoided.
- This product is intended for use on one person. It may not be reused on another person.

Specific Safety Instructions

- Orthosis cannot be thermoformed.
- Do not bring orthosis into contact with acidic and fatty substances or solutions.
- Avoid local pressure points and possible concurrent disruption of blood flow by not fastening the orthosis too tighly.
- The orthosis can only be used for a limited range of athletic activities.
 If the patient wants to perform any athletic activities, the orthopedic technician should inform her or him

- about the loading capacity and resilience of the orthosis.
- Movements during which the orthosis gets bent excessively (e.g. initial contact of foot with the whole forefoot instead of the tip of the toes while climbing/descending stairs, or squatting) should be avoided.
- Maximum body weight is 100 kg.

REF/Sizes

Return-to-Walk P, Medium Stiffness				
REF Left	REF Right	Size	Shoe Size	Height of Orthosis
22.390.302.00	22.390.402.00	S	35-38	35 cm
22.390.303.00	22.390.403.00	М	38-41	36 cm
22.390.304.00	22.390.404.00	L	41-44	37 cm
22.390.305.00	22.390.405.00	XL	44-47	37 cm

Return-to-Walk P/Plus, High Stiffness				
REF Left	REF Right	Size	Shoe Size	Height of Orthosis
22.390.502.00	22.390.602.00	S	35-38	35 cm
22.390.503.00	22.390.603.00	М	38-41	36 cm
22.390.504.00	22.390.604.00	L	41-44	37 cm
22.390.505.00	22.390.605.00	XL	44-47	37 cm

Scope of Delivery

Carbon orthosis with calf pad and sewn in hook-and-loop fastener, Included are a pad (optional use) to protect the medial ankle. an insole with two adhesive tapes and a user manual with all relevant information.

the medial strut is positioned behind the medial ankle and does not touch it. Carefully trim the forefoot area as necessary. Make sure not to sand in the area of the carbon fibers, (fig. 2).

Application/Adjustment

- **1.** Determine the shoe size of the patient and select the size of the orthosis accordingly.
- 2. For best results the patient should be wearing a sturdy outdoor shoe with a strong heel counter. Heel height should he about 1 cm
- **3.** If the insole of the patient's shoe can be removed, use it as a pattern to transfer the right size to the foot plate (fig. 1).

4. Only adjust the width of the foot plate on the lateral side, do not adjust more than required. Adjust the length of the foot plate by shortening the heel area, so

5. The textile pad should lie against the top calf part. Cut the calf strap with the hook-and-loop fastener to the right length, so the end of the hookand-loop fastener is positioned in the extension of the strut (fig. 3a/3b).

Comfort/Hygiene

Do not wear the orthosis directly on the skin, but instead over thin clothing (e.g. a sock of appropriate length).

Material

The orthosis is made from carbon fiber composite material. The textile component is made of velours, soft hookand-loop material and spacer fabric. The product is latex-free. The insole is made from Evazote foam material.

Produkt Care/ Cleaning/Maintenance

• Textile Component:

Hand wash only.

Do not bleach.

Do not tumble dry.

Do not iron.

Do not dry-clean.

• Carbon Fiber Composite Component: Clean with a wet cloth.

Storage/Disposal

Store the product in a dry place. Protect from dust, dirt and moisture. This product does not contain any environmentally hazardous substances. It can be disposed of with household waste.

CE-Conformity

This medical device complies with the requirements of the MDR, Regulation (EU) 2017/745.

Legend to Symbols

Manufacturer

Date of Manufacture

Distributor

Importer

EC Authorized Representative

Lot Number

Part Number

Serial Number

UDI Data Matrix

Follow User Manual

Caution

CE-Marking According to European Union Council Regulation 2017/745 Concerning Medical Devices (MDR)

Single Patient, for Multiple Use

Store in a Dry Place, Protect from Moisture

Protect from Sunlight

Temperature Limit

Humidity Limit

Atmospheric Pressure Limit, for Transport and Storage

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Wilhelm Julius Teufel GmbH Robert-Bosch-Straße 15 73117 Wangen Deutschland/Germany