
TKS Knieorthese

TKS Knee Orthosis

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User Manual

◆ Inhalt

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Content

Description	EN-2
Indication	EN-2
Contraindication	EN-2
Function	EN-2
Reimbursement Code	EN-2
Medical Device	EN-2
Safety	EN-3
REF/Sizes	EN-3
Taking Measurements	EN-4
Scope of Delivery	EN-4
Application/Adjustment	EN-4
Comfort/Hygiene	EN-4
Material	EN-4
Product Care/Cleaning/Maintenance	EN-4
Storage/Disposal	EN-4
CE-Conformity	EN-5
Legend to Symbols	EN-5

Description

The TKS knee orthosis helps control structural instabilities of the knee. The construction of the brace is based on the three-point-principle, which stabilizes the knee joint, allows for guided flexion and at the same time prevents hyperextension. The orthosis consists of two lateral aluminum braces which are linked by two aluminum bands. The padded lining of the orthosis offers a high level of wearing comfort. On the front of the leg, an aluminum bands runs across the thigh above the knee, on the back of the leg, the second aluminum bands runs across the calf below the knee. The upper fastening is positioned at the back of the thigh. It is closed with a strong textile strap. The bottom fastening is

positioned on the front of the lower leg. It features an anatomically shaped, cushioned pad which helps avoid too much pressure on the tibia. The knee joint's range of motion can be manipulated by adjusting the hinged connections and by adjusting the fastenings.

Indication

Genu recurvatum

Contraindication

- ACL instabilities
- Flexion contractions

Function

- Stabilizing the knee joint according to the three-point-principle
- Effective, adjustable control of genu recurvatum

Reimbursement Code

DE: HMV-Nr.: 23.04.07.0003

◆ Medical Device

The TKS Knee Orthosis is a medical device.

◆ **A** Safety

- Before this device is used, the orthosis has to be individually adapted to the patient by an experienced orthopedic technician and the patient has to be trained in the orthosis' use according to the doctor's instructions.
- Incorrect selection, adaptation, application, and/or failure to monitor the correct fit of this product may cause health damage.
- Check orthosis for damage (such as loose connections, fasteners, belts, etc.) each time it is put on.
- Do not remove the textile label from textile products.
- In case you suspect any performance malfunctions (e.g. an apparent reduction or lack of effect), immediately call on your doctor or orthopedic technician.
- All serious incidents that occur in connection with this product have to be reported to the Wilhelm Julius Teufel GmbH and the competent authority of the member state in which you are established. A "serious incident" is any incident that caused, may have caused, or may cause,

- directly or indirectly, one of the following consequences:
- a) death of a patient, user, or other persons
- b) temporary or permanent health deterioration of a patient, user or other persons
- c) a serious public health threat.
- Never wear the orthosis on broken skin.
- The orthosis must be checked for proper function, correct fit and potential damage by a doctor or orthopedic technician at regular intervals. The timing of these intervals has to be defined individually for each patient.
- Should the product have been subject to disproportionate stress, it must be checked for potential damage by an orthopedic technician before further use.
- Harmful pressure caused by the orthosis may cause skin injuries or disrupt blood flow and must be avoided.
- This product is intended for use on one person. It may not be reused on another person.

REF/Sizes

REF Left Knee	REF Right Knee	Size	Clear Width*	Thigh Circumference
21.277.103.00	21.277.203.00	3	90 mm	up to 40 cm
21.277.104.00	21.277.204.00	4	100 mm	41-45 cm
21.277.105.00	21.277.205.00	5	110 mm	46-50 cm
21.277.106.00	21.277.206.00	6	120 mm	51-55 cm
21.277.107.00	21.277.207.00	7	130 mm	56 cm and up

^{*} The lateral aluminum braces are 240 mm long for all sizes.

Taking Measurements

Select the required Size based on the thigh circumference 15 cm above the Patella. Alternatively, measure the clear width at the height of the knee joint gap.

Scope of Delivery

REF	Produkt
21.277.x0x.00	TKS Knee-Orthosis
90.395.018.01	TKS_UserManual_ DE_EN

◆ Application/Adjustment

- 1. Open all hook and loop fastenings.
- 2. Step through the orthosis with the affected leg.
- 3. Position the thigh band 1-2 cm above patella. Make sure that the lateral uprights are parallel to the leg.
- 4. Close the strap on the back of the thigh. Once the length of this strap has been adjusted, the strap should only be opened and closed using the lateral snap button. The hook-and-loop fastening should not be opened anymore.
- 5. Place the anatomically shaped pad onto the front of the lower leg. The cut-out should be positioned on the tibial edge. Close the strap.
- 6. The two-hole-nuts fastening the aluminum straps to the aluminum braces have to be secured with threadlocker after fitting the orthosis.

◆ Comfort/Hygiene

Do not wear the orthosis directly on skin, but instead over thin clothing (e. g. a pair of leggings).

Material

Aluminium Polyethylene Foam Textile

Product Care/ Cleaning/Maintenance

Use a damp cloth to clean the TKS Knee Orthosis

Hand wash.

Do not bleach.

Do not tumble dry.

Do not iron.

Do not dry clean.

◆ Storage/Disposal

The orthosis does not contain any environmentally hazardous substances. It can be disposed of with household refuse.

CE-Conformity

This medical device complies with the requirements of the MDR, Regulation (EU) 2017/745.

Legend to Symbols

Manufacturer

Date of Manufacture

Distributor

Importer

EC Authorized Representative

Lot Number

Part Number

Serial Number

UDI Data Matrix

Follow

User Manual

Caution

CE-Marking According to European Union Council Regulation 2017/745 Concerning Medical Devices (MDR)

Single Patient, for Multiple Use

Store in a Dry Place, Protect from Moisture

Protect from Sunlight

Temperature Limit

Humidity Limit

Atmospheric Pressure Limit, for Transport and Storage

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Wilhelm Julius Teufel GmbH Robert-Bosch-Straße 15 73117 Wangen Deutschland/Germany