CoxaFlex®

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User Manual
Modo de empleo
Manuale d'instruzioni

CE

www.wjt-ortho.com

WJT
ORTHO

◆ Inhalt

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◆ Contents

Description	EN-2
Indication	EN-3
Contraindication	EN-3
Function	EN-3
Reimbursement Code	EN-3
Medical Device	EN-3
Safety	EN-3
Intended Product Use	EN-4
REF/Sizes	EN-5
Scope of Delivery	EN-5
Taking Measurements	EN-6
Application/Adjustment	EN-6
Comfort/Hygiene	EN-8
Material	EN-8
Product Care/Cleaning/Maintenance	EN-8
Storage/Disposal	EN-8
CE-Conformity	EN-8
Legend to Symbols	EN-8

position). Evasive movements leading to a wider flexion and rotating motions, however, are deliberately promoted. Therefore, the flexion abduction splint CoxaFlex allows for a precisely controllable, unconstrained adjustment of the lower extremities in flexion and abduction. It serves the centering of the femoral head with an instable or dysplastic hip joint structure (i.e. developmental hip disorder) under circumstances as natural as possible. CoxaFlex helps maintain the foetal "seated squat position" recommended for newborns and infants and also complies with the natural squatting position of children huddling against their mother's body.

Where applied correctly, no unwanted side effects have been reported for CoxaFlex in more than 25 years.

The CoxaFlex flexion abduction splint is not capable of reducing a luxated hip, but of centering and stabilizing the joint in a physiological position. To prevent relaxation, special attention should hence be paid to the hip joint's correct and reduced position when applying the orthosis.

In case of dysplasia (without the risk of luxation) the orthosis can shortly be taken off during the day (e.g. for hygienic care). This, however, should only be done with the explicit permission of the attending physician.

If worn over an opening romper, the orthosis does not have to be taken off when changing diapers. The splint's convex shaping in the belly area allows for the diaper to be easily removed.

◆ ⚠ Safety

◆ Indication

- Securing hip joints that are prone to dislocate
- Hip dysplasia
- Set hip luxations

◆ Contraindication

Non-reduced hip luxation.

◆ Function

- Alignment of the hip joints
- Adjustment of the hip joints in abduction and flexion

◆ Reimbursement Code

DE HMV-Nr.: 23.05.02.0012

AT Tarif Nr.: 35033 019

◆ Medical Device

The CoxaFlex is a medical device.

- Before this device is used, the orthosis has to be individually adapted to the patient by an experienced orthopedic technician and the patient has to be trained in the orthosis' use according to the doctor's instructions.
- Incorrect selection, adaptation, application, and/or failure to monitor the correct fit of this product may cause health damage.
- Check orthosis for damage (such as loose connections, fasteners, belts, etc.) each time it is put on.
- Do not remove the textile label from textile products.
- In case you suspect any performance malfunctions (e.g. an apparent reduction or lack of effect), immediately call on your doctor or orthopedic technician.
- All serious incidents that occur in connection with this product have to be reported to the Wilhelm Julius Teufel GmbH and the competent authority of the member state in which you are established. A „serious incident“ is any incident that caused, may have caused, or may cause, directly or indirectly, one of the following consequences:
 - a) the death of a patient, user, or other persons
 - b) temporary or permanent health deterioration of a patient, user or other persons
 - c) a serious public health threat.
- Never wear the orthosis on broken skin.
- The orthosis must be checked for proper function, correct fit and potential damage by a doctor or orthopedic technician at regular intervals. The timing

of these intervals has to be defined individually for each patient.

- Should the product have been subject to disproportionate stress, it must be checked for potential damage by an orthopedic technician before further use.
- Harmful pressure caused by the orthosis may cause skin injuries or disrupt blood flow and must be avoided.
-  This product is intended for use on one person. It may not be reused on another person.

◆ Intended Product Use

CoxaFlex was developed in order to treat hip dysplasia and luxation more effectively and more safely. It helps control the growth of the cartilaginous and bony articulation partners, and specifically supports the ideal shaping of the hip joint.

CoxaFlex resulted from the close cooperation of physicians and technicians and will foster the subsequent development of the hip joint when applied correctly.

As with many disorders, the prospects of a complete recovery of the dysplastic hip will improve if the defect is detected and treated early on.

When treated within the first days and weeks of life, regular dysplasia can usually be healed within a few weeks. Dysplasia stemming from a preceding luxation, however, usually takes several months of orthotic treatment to heal.

Constant and regular control examinations (at first every 2 to 3 days) are as much part of a good treatment plan as a professionally fitted and correctly applied orthosis. Therefore, the close cooperation of physicians, technicians and parents is indispensable for the patient's healthy future.



Newborns and infants show a high degree of strong, reflexive, kicking leg movements. Thus, children subconsciously reject any form of restriction to their mobility, and might show some resistance against the splint until fully

accustomed to the new orthosis.

Should a patient seem to cry out of pain for an unusually long period of time, remove the orthosis and see a physician

so the proper fit of the splint and the stability of the joints can be checked.

Children wearing the splint should preferably sleep on their backs.


◆ REF/Sizes

REF	Size	Min. Brace Spread/cm	Max. Brace Spread/cm	Weight/g
23.300.000.00	XS	8.0	13.5	148
23.300.001.00	S	9.5	16	153
23.300.002.00	M	11	17.5	172
23.300.003.00	L	12.5	19	192
23.300.004.00	XL	17.5	24	200

◆ Scope of Delivery

REF	Product
23.300.00X.00	CoxaFlex Flexion Abduction Splint
90.395.070.01	CoxaFlex_UserManual_DE_EN_ES_IT

◆ Taking Measurements

◆ Application/Adjustment

1. Put the orthosis on a flat surface with the blue neck part facing you. Open the white buckle in the chest area and open the two front buckles on the leg braces (fig. 1).

2. Completely unfold the orthosis in front of you (fig. 2).

3. Lay the patient onto the belt system and put the straps over the infant's shoulders. Place the child's thighs into the leg braces (fig. 3).

4. First, close the belt buckles on the leg braces, then the white buckle on the chest. The flexion abduction splint is usually applied over the infant's romper. (Exception: With unstable hips the splint should be worn under the romper, with tubular cotton stockinette placed underneath the orthosis to protect the skin.) During the fitting, have a second person hold the child's legs so the hip joints are bent at 100° and are abducted at an


angle of 50°. Now adjust the distance between the knees (fig. 4).

5. Position the closure of the chest belt as low as possible (fig. 5).

6. Adjusting and securing the belts: Pull the belts out of the three-bar-slides and adjust them to the desired length before putting them back into the slides. The belts should be suspended in such a way that a flexion of 90° to 100° can be ensured even during strong pedalling motions (fig. 6).

⚠ Please Note: The knurled screw generally allows for a completely variable adjustment of the degree of abduction without using any tools. To ultimately secure the desired adjustment, the knurled screw can be fastened using a coin. The widest brace spread can be achieved by extending the perforated rails to the maximum. Always verify that the buckles on the leg braces are securely locked (fig. 7).

◆ Comfort/Hygiene

Do not wear the orthosis directly on the skin, but instead over clothing.

◆ Material

- Stainless steel
- Polyester
- Subortholen

◆ Product Care/ Cleaning/Maintenance

- Neck and brace cushions are removable and should be washed with a detergent for delicate fabrics at 30°C. A spare neck cushion is included free of charge with every CoxaFlex.
- Neck part, belts and spreading-unit should only be wiped down with a damp cloth.

Wash on delicate cycle at 30°C.

Do not bleach.

Do not tumble dry.

Do not iron.

◆ Storage/Disposal

Store the product in its original packaging in a dry place. Protect from dust, dirt and moisture. This product does not contain any environmentally hazardous substances. It can be disposed of with household refuse.

◆ CE-Conformity

This medical device complies with the requirements of the MDR, Regulation (EU) 2017/745.

◆ Legend to Symbols

Manufacturer

Date of Manufacture

Distributor

Importer

EC Authorized Representative

Lot Number

Part Number

Serial Number

UDI Data Matrix

Follow User Manual

◆ Notes

Caution

CE-Marking According to European Union Council Regulation 2017/745 Concerning Medical Devices (MDR)

Single Patient,
for Multiple Use

Store in a Dry Place,
Protect from Moisture

Protect from Sunlight

Temperature Limit

Humidity Limit

Atmospheric Pressure Limit,
for Transport and Storage

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