
DynaCox® evolution Hüftorthese

Hip Orthosis

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User Manual

Inhalt

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Contents

Description	EN-2
Indication	EN-2
Contraindication	EN-2
Function	EN-2
Reimbursement Code	EN-2
Medical Device	EN-3
Safety	EN-3
Intended Product Use	EN-4
REF / Sizes	EN-4
Taking Measurements	EN-5
Scope of Delivery	EN-7
Accessories	EN-7
Application / Adjustment	EN-8
Comfort / Hygiene	EN-20
Material	EN-20
Product Care / Cleaning / Maintenance	EN-20
Storage / Disposal	EN-20
CE-Conformity	EN-20
Legend to Symbols	EN-21

Description

The DynaCox evolution hip orthosis provides extensive luxation protection for the hip joint. This is achieved through the stabilizing effect of an orthosis working according to the three-point-principle, with joint control, combined with a functional bandage.

The elastic, fully-fashioned knit fabric ensures an even compression and optimal wearing comfort.

You can choose from two joint versions:

 On the ROM-Flex joint, flexion and extension can be adjusted in increments of 15°. • The TRI-Flex joint is limited to a maximum flexion of 70°. Extension is not limited.

The modular design of the orthosis allows for quick application and full adaptability.

The DynaCox evolution hip orthosis is approved for clinical and domestic use.

Indication

- After repositioning of luxated hip joints
- Instability after Total Hip Replacement (THR)
- Muscular imbalance

Contraindication

- Non-compliant patients
- Severe skin alterations

Function

- Securing surgical results
- Luxation prophylaxis
- Early functional motion therapy

Reimbursement Code

(DynaCox evolution with ROM-Flex Joint)

DE HMV-No.: 23.05.01.0002 **AT** Tarif-No.: 35333 020

Medical Device

The DvnaCox evolution is a medical device.

Safety

- Before this device is used, the orthosis has to be individually adapted to the patient by an experienced orthopedic technician and the patient has to be trained in the orthosis' use according to the doctor's instructions.
- Incorrect selection, adaptation, application, and/or failure to monitor the correct fit of this product may cause health damage.
- Check orthosis for damage (such as loose connections, fasteners, belts, etc.) each time it is put on.
- Do not remove the textile label from textile products.
- In case you suspect any performance malfunctions (e.g. an apparent reduction or lack of effect), immediately call on your doctor or orthopedic technician.
- All serious incidents that occur in connection with this product have to be reported to the Wilhelm Julius Teufel GmbH and the competent authority of the member state in which you are established. A "serious incident" is any incident that caused, may have caused, or may cause, directly or indirectly, one of the following consequences: a) death of a patient, user, or other persons

- b) temporary or permanent health deterioration of a patient, user or other persons
- c) a serious public health threat.
- Never wear the orthosis on broken skin.
- The orthosis must be checked for proper. function, correct fit and potential damage by a doctor or orthopedic technician at regular intervals. The timing of these intervals has to be defined individually for each patient.
- Should the product have been subject to disproportionate stress, it must be checked for potential damage by an orthopedic technician before further use.
- Harmful pressure caused by the orthosis may cause skin injuries or disrupt blood flow and must be avoided.
- This product is intended for use on one person. It may not be reused on another person.

The DynaCox evolution hip orthosis is designed to be adjusted for either the right or the left hip joint, but not for bilateral use. Due to the required individual adaptation the orthosis is not suitable for reuse by other persons!

◆ Intended Product Use

This product is intended to be used on one patient only.

Compatibility

The DynaCox evolution may not be combined with products other than those mentioned in this user manual.

Weaning Off

Following the course of the therapy, the rigid (orthotic) components of the DynaCox evolution can be separated from the bandage.

◆ REF / Sizes

DynaCox® evolution with TRI-Flex Joint, including bandage					
REF Left	REF Right	Size	Circumference/cm		
			Hip	Thigh*	
23.600.101.00	23.600.201.00	S	70 - 85	40 - 52	
23.600.102.00	23.600.202.00	Μ	85 - 100	48 - 60	
23.600.103.00	23.600.203.00	L	100 - 120	56 - 68	
23.600.104.00	23.600.204.00	XL	115 - 135	62 - 76	
23.600.105.00	23.600.205.00	XXL	130 - 150	65 - 80	

^{*}Refer to thigh circumference in case the patient's hip circumference is close to the next bigger or next smaller size.

DynaCox® evolution with ROM-Flex Joint, including bandage					
REF Left	REF Right	Size	Circumference/cm		
			Hip	Thigh*	
23.600.101.01	23.600.201.01	S	70 - 85	40 - 52	
23.600.102.01	23.600.202.01	Μ	85 - 100	48 - 60	
23.600.103.01	23.600.203.01	L	100 - 120	56 - 68	
23.600.104.01	23.600.204.01	XL	115 - 135	62 - 76	
23.600.105.01	23.600.205.01	XXL	130 - 150	65 - 80	

^{*}Refer to thigh circumference in case the patient's hip circumference is close to the next bigger or next smaller size.

◆ Taking Measurements

- 1. The starting point for measuring in the pelvis area is the greater trochanter.
- 2. Measure the hip circumference (about 6 cm above the greater trochanter, fig. 1).

3. Write down the measurement.

4. Please refer to the table below as a guide when adjusting the width of the pelvic brace.

Size Adjustment Pelvic Brace					
Hip Circumference (cm), Orthosis with Bandage	Hole on Posterior Belt, Counted from the Middle of the Belt	Hole on Anterior Belt, Counted from the Middle of the Belt (Buckle)			
76	1	1			
78	1	1			
81	2	1-2			
84	2	1-2			
87	3	1-2			
90	3	1-2			
93	4	2			
96	4	2			
99	5	2-3			
102	5	3			
105	6	3-4			
108	6	4			
111	7	4-5			
114	7	5			
117	8	5-6			
120	8	6			
123	9	6-7			
126	9	7			
129	10	7-8			
132	10	8			
135	11	8-9			
138	11	9			
141	12	9-10			
144	12	10			
147	13	10-11			
150	13	10-12			

Please note: The numbers in the 2nd and 3rd column denote the hole in the belt into which the outer peg of the pelvic panel should be locked to achieve a good fit.

Scope of Delivery

Bandage, Joint, Solid Components, Spare Parts Kit, Tool Kit, User Manual, Order Form

Accessories

Bandage for DynaCox® evolution, replacement bandage

DynaCox® evolution	DynaCox® evolution Bandage, Right	Size	Circumference/cm	
Bandage, Left			Hip	Thigh*
23.601.101.00	23.601.201.00	S	70 - 85	40 - 52
23.601.102.00	23.601.202.00	М	85 - 100	48 - 60
23.601.103.00	23.601.203.00	L	100 - 120	56 - 68
23.601.104.00	23.601.204.00	XL	115 - 135	62 - 76
23.601.105.00	23.601.205.00	XXL	130 - 150	65 - 80

We recommend the use of a spare bandage (REF 23.601.XXX.XX) and a spare set of pads (REF 23.604.020.00) when the DynaCox evolution has to be worn for an extended period of time.

^{*}Refer to thigh circumference in case the patient's hip circumference is close to the next bigger or next smaller size.

REF	Product
23.603.000.00	DynaCox® evolution ROM-Flex-Joint
23.603.100.00	DynaCox® evolution TRI-Flex-Joint Left
23.603.200.00	DynaCox® evolution TRI-Flex-Joint Right

REF	Product	Scope of Delivery
23.604.010.00	DynaCox evo Solid Components f.Consignm.	DynaCox evolution Rigid Components, DynaCox evolution Tool Kit, DynaCox evo Spare Parts Kit
23.604.020.00	DynaCox evo Pad Set	Condyle Pad Left/Right, Condyle Clasp Pad Air, Condyle Clasp Pad Foam, Pelvic Brace Pad Left, Pelvic Brace Pad Right, Pelotte Padding, Abdominal Textile Strap, Condyle Clasp Textile Strap
23.604.030.00	DynaCox evo Strap Set	Pelvic Strap Plastic, Abdominal Strap Plastic, Buckle f. 35mm-Straps Black Lateral
23.604.040.00	DynaCox evo Lateral Bar, short, with pad	Short Bar, Condylar Pad (For patients shorter than 1.55 m)
23.604.041.00	DynaCox evo Lateral Bar, standard, with pad	Standard Bar, Condylar Pad
23.604.050.00	DynaCox evo Panel Padding, with hook and loop fastener	Panel Padding (For exceptionally slender patients)
23.604.060.00	DynaCox evo Spare Part Set	Cylinder Head Screw, Button Head Screw, Spring Washer, Washer
60.397.001.00	Heavy Duty Scissors	1 pair of heavy duty scissors, serrated.
23.600.990.00	DynaCox evolution Tool Kit	Multi-Purpose Scissors, Hex-Wrench 3x75mm

Application / Adjustment

Please note: The DynaCox evolution can be applied with the patient either standing or lying. Wherever necessary, the following instructions will point out the differences between these two types of application.

Pre-Adjusting the **Orthosis Component**

To keep the application time as short as possible for the patient, please make all adjustments on the orthosis component without the bandage. Only connect the orthosis component to the bandage

when all of the following pre-adjustments have been completed.

Pre-Adjusting the Pelvic Brace

1. The pelvic brace consists of four components (see fig. 2):

• an anterior belt (recognizable by the central buckle)

EN-8

- a posterior belt (recognizable by a sticker saving "UP/INSIDE")
- two pelvic panels

The width of the hip brace can be adjusted by changing the length of the anterior and posterior belt. Use the hip measurement as an orientation when initially adjusting the belts.

2. Remove the pads from the pelvic panels, and adjust the posterior belt so that the lateral Click & Lock closures on the pelvic panels are positioned on the lateral femoral center line (see fig. 3).

- **3.** Detach the posterior belt from the two pegs on the pelvic panel.
- **4.** Pull or push the posterior belt to the desired length.
- **5.** Snap the pegs on the pelvic panel back into the holes on the belt to secure the adjustment. Make sure to first snap the inner, then the outer peg nto place (fig. 4).

- **6.** Repeat the process on the other end of the posterior belt.
- 7. Now adjust the anterior belt. Detach the belt on one side from the two pegs on the pelvic panel.
- 8. Pull or push the anterior belt to the desired length.
- **9.** Snap the pegs on the pelvic panel back into the holes on the helt to secure the adjustment. Please make sure to first snap the inner, then the outer peg into place.
- **10.**Repeat the process on the other end of the anterior belt.
- **11.** If the patient is very slender it can be helpful to shorten the posterior and anterior belt on the inside of the pelvic panels after attaching them to the panels. To do so, cut the belts with a suitable pair of scissors (REF 60.397.001.00) (fig. 5).

Please make very sure that the pelvic brace has been adjusted properly before shortening the belts.

Adjusting the ROM-Flex Joint

Pull the black joint cover back, far enough so you are able to remove the screws (fig. 6).

Please note: Always open only one of the stop screws at a time. If both screws are opened at the same time, the threaded ring on the rear of the joint will come undone. However, the rear ring can be screwed back on, without any loss of functionality to the joint.

1. Upon delivery, the ROM-Flex joint is set to 60° flexion and extension (fig. 7).

- **2.** To change the range of motion, completely unscrew one of the stop screws.
- **3.** The current angle value can be seen on the engraved scale on the joint.

- **4.** Place the stop screw in one of the threads, so that the desired flexion or extension can be achieved. Now tighten the screw again.
- **5.** Flexion and extension can be adjusted in increments of 15°.
- **6.** Unscrew the second stop screw and repeat the process.
- Tighten this second stop screw again in order to finish the adjustment of the joint.
- **8.** Now push the black joint cover back over the ring of the joint, covering the screws completely.

TRI-Flex-Joint

The TRI-Flex joint is available as a left or right joint. The joints are marked with the letter "L" or the letter "R". They may only be installed on the indicated side, since otherwise flexion will not be limited. The TRI-Flex joint is limited to a maximum flexion of 70°. Extension is not limited.

Adjusting the Length of the Lateral Bars

Due to the Click & Lock closure it is possible to adjust the length of the joint in relation to the pelvic panel.

- To adjust the length of the joint, insert the bar of the joint (TRI-Flex or ROM-Flex) into the open Click & Lock closure.
- **2.** It is absolutely mandatory not to exceed the maximum extension

- length. Otherwise, a secure connection between the bar and the pelvic panel cannot be guaranteed.
- **3.** The bar is secured in its position by the pins on the Click & Lock closure. which lock into the holes in the bar. In case the Click & Lock closure cannot be locked, push or pull the bar up or down a little within the limits of the adjustment area, so the pins can lock into the holes on the har.
- **4.** Double Pin: The double pin version can be identified by the two metal pins, which have been embedded directly into the lever of the Click&Lock closures (fig. 8).

Please note: The maximum permitted extension length of the bar. There is an engraved line on the bar that indicates the maximum permitted extension. The bar must not be pulled out of the Click & Lock closure beyond this line. (Fig. 9: maximum extension length.)

- 5. Hold the orthosis alongside the patient's affected leg for orientation.
- **6.** Now insert the upper lateral bar into the Click & Lock closure on the pelvic panel.
- 7. Adjust the length of the bar so that the upper edge of the femoral panel sits directly under the greater trochanter. At the same time, the pivot point of the orthosis' joint should be positioned congruently with the anatomical pivot point of the hip joint (fig. 10).

8. To secure the upper bar's position, only push the Click & Lock closure in until it touches the pelvic panel, but do not push the lever all the way into the cut-out on the panel. Do not lock the closure completely yet, so it remains possible to make adjustments later on, if necessary.

Adjusting the Condyle Clasp

1. The DynaCox evolution comes with two different pads for the condylar clasp, i.e. a thicker and a thinner one (fig. 11). The thicker pad is cushioned with air pockets. Both pads feature markings for potential cutting lines along the edges. If the patient's leg circumference is very small, it may be necessary to shorten the pads on both sides, in order to be able to adjust the condylar pads as necessary. We recommend that the thinner pad be used for obese patients.

2. Upon delivery, the lateral bar of the condylar clasp is already installed on the designated side. Please note that if the products were ordered for a consignment stock the lateral side may have to be transferred to the other side. If you want to switch the side on which the bar is attached, remove the two screws on the side of the condylar clasp, where the lateral bar is attached. Now pull the condylar clasp off the condylar pad (fig. 12), turn it by 180°, and reattach it on the condylar clasp (fig. 13). Then screw the condylar clasp and the condylar pad back together.

3. Now adjust the position of the condylar pads. Do so by pushing the pads along the slide track. This way, the pads can be moved into a more anterior or posterior position (fig. 14). Secure the position by carefully tightening the two hex screws (1Nm).

4. In case your patient's leg circumference is very large, it may be necessary to adjust the angle of the condylar clasp to the angle of the lateral bar in order to achieve a good fit. To do this, loosen the two hexagon socket screws connecting the lateral bar with the condylar clasp. Take care to loosen, but not to remove the screws. It is now possible to tilt the condylar clasp by

7° (fig. 15/16). After the position of the condylar clasp has been adjusted, tighten the screws again (1Nm).

The condylar pads should rest against the leg on the medial and lateral side, above the condyles (fig. 17).

5. Insert the lower lateral bar into the Click & Lock closure on the femoral panel (fig. 18). At the most, the lateral bar can be pushed into the femoral panel until it reaches the top stop. A short lateral bar, which can be used if the patient's femur is exceptionally short, is optionally available (REF 23.604.040.00). The lateral bars must not be shortened under any circumstances!

- **6.** To secure the bar's position, only push the Click & Lock closure's lever in until it iust touches the femoral panel, but do not push the lever all the way into the cut-out on the panel.
- 7. Do not lock the Click & Lock closure completely vet, so it remains possible to make adjustments later on, if necessarv. In case the closure can later not be locked completely, try to push or pull the bar up or down a little while the Click & Lock closure is already pushed in, so the pins can audibly lock into a hole.
- 8. Put the pads back into the pelvic panels.

Attaching the Orthosis Component to the Bandage

- 1. Open the central buckle on the anterior helt
- **2.** Attach the orthosis component to the bandage by pushing the pelvic brace through the loop on the back of the bandage (fig. 19).

3. Now insert the femoral panel into the pocket on the bandage (fig. 20).

4. Close the zipper of the bandage over the lower lateral bar (fig. 21).

5. Close the hook-and-loop fastener on the upper end of the pocket on the bandage.

Technician's Tipp:

In case there is a high degree of atrophy on the patient's leg, doubling the padding on the femoral panel may be helpful to achieve a proper fit. An additional pad can be ordered with REF 23.604.050.00 (fig. 32).

• Fitting the Orthosis on the Patient

Fitting the Orthosis on a Standing Patient

Please note: The patient should be standing with slightly abducted legs to allow for easier adjustment of the DynaCox evolution.

- 1. When putting the brace on the patient, the upper edge of the femoral panel should be placed directly below the greater trochanter, and the lateral bars should be positioned on the lateral femoral center line.
- **2.** Close the bandage's pelvic hook-and-loop closure (fig. 22).

3. Now put the condylar clasp around the patient's leg, and check the position of the condylar pads. The position of the condylar pads can be adjusted step by step, by pushing the pads further to the front or back of the condylar clasp (fig. 23).

4. Secure the position of the condylar pads (1Nm), by lightly tightening the hex head screws (fig. 24).

6. Close the condylar clasp by feeding the hook-and-loop strap into the strap guide and closing it (fig. 27).

5. Now close the pelvic brace of the orthosis component (fig. 25).

7. Close the thigh section of the bandage (fig. 28).

If necessary, adjust the length of the anterior and posterior belt again. Minimal adjustments can also be made with the straps to the right and the left of the central buckle on the anterior belt (fig. 26). If necessary, shorten the straps.

The loop strap of the hook-and-loop closure on the bandage can be shortened with a suitable pair of scissors (REF 60.397.001.00). Make sure to cut between the loop fabric chevrons when shortening the strap.

- **8.** Re-check the correct fit of the orthosis, and make some more adjustments, if necessary.
- **9.** Once a correct fit has been achieved, secure the Click & Lock closures on the femoral panel and the hip brace.

Please note: Only close the Click & Lock closures by pushing against the latch with your thumbs while <u>simultaneously</u> creating counter pressure with your fingers. This will ensure that neither the leg nor the hip of the patient are subjected to any possibly harmful pressure (fig. 29).

- 10. It may be necessary to open the central buckle of the pelvic brace again so you can easily push the latches of the Click & Lock closures into position. Afterwards, the central buckle can be closed again.
- 11. Please consider the following when closing the Click & Lock closure on the panel: For all length settings (minimum, maximum and in between) each of the two pins must lock into one of the holes in the bar. Only then the Click & Lock closures will sit flush against the orthosis.

- **12.** In case the Click & Lock closure can later not be locked completely, push or pull the bar up or down a little, until both pins have audibly locked into the next hole.
- **13.** After pushing the latches flush against the orthosis, any unintended opening of the DynaCox evolution will not be possible anymore. To open the Click & Lock closures, insert a hex key into the small opening on the back of the closure (fig. 30), and push the latch open (fig. 31). This opening can be found on both the closures on the femoral panel and the pelvic panels.

Technician's Tipp:

In case there is a high degree of atrophy on the patient's leg, doubling the padding on the femoral panel may be helpful to achieve a proper fit (fig.32). An additional pad can be ordered with RFF 23.604.050.00.

Fitting the Orthosis on a Lying Patient

Please note: Especially when fitting the orthosis on a lying patient, accurately preadjusting the DynaCox evolution is very important in order to make the process as simple as possible for both patient and orthopedic technician.

- **1.** During the application, the patient should be lying on her/his back.
- **2.** Pass the pelvic brace together with the bandage under the patient's lumbar region (fig. 33).

- 3. Position the DynaCox evolution in a way that places the upper edge of the femoral panel directly under the greater trochanter, so that the lateral bars are positioned on the lateral femoral. center line.
- **4.** First close the bandage's pelvic hook and loop closure (fig. 34).

5. Now close the pelvic brace of the orthosis component with the central buckle (fig. 35).

Minor adjustments can also be made with the straps to the right and the left of the central buckle on the anterior belt (fig. 36).

6. Put the condylar clasp around the patient's leg, and check the position of the condylar pads. The position of the condylar pads can be adjusted step by step, by pushing the pads further to the front or to the back of the condylar clasp (fig. 37).

7. Secure the position of the condylar pads by carefully tightening the hex head screws (1Nm) (fig. 38).

8. Close the condylar clasp by feeding the hook-and-loop strap into the strap guide and closing it (fig. 39).

9. Close the thigh section of the bandage (fig. 40).

The loop strap of the hook-and-loop closure on the bandage can be shortened with a suitable pair of scissors (REF 60.397.001.00). Make sure to cut between the loop fabric chevrons when shortening the strap (fig. 41).

- **10.** Re-check the correct fit of the orthosis. Make more adjustments if necessary.
- **11.** Once a correct fit has been achieved, secure the Click & Lock closures on the femoral panel and the hip brace.

Please note: Only close the Click & Lock closures by pushing against the latch with vour thumbs while simultaneously creating counter pressure with your fingers (fig. 42). This will ensure that neither the leg nor the hip of the patient are subjected to any possibly harmful pressure.

- **12.** It may be necessary to open the central buckle of the pelvic brace again in order to be able to easily push the latches of the Click & Lock closures into position. Afterwards, the central buckle can be closed again.
- 13. Please consider the following when closing the Click & Lock closure on the panel: For all length settings (minimum, maximum and in between) each of the two pins must lock into one of the holes in the bar. Only then the Click & Lock closures will sit flush against the orthosis. In case the Click & Lock closure can later not be locked. completely, push or pull the bar up or down a little, until both pins have audibly locked into the next hole.

14. After pushing the latches flush against the orthosis, any unintended opening of the DynaCox evolution will not be possible anymore. To open the Click & Lock closures, insert a hex key into the opening on the back of the closure (fig. 43), and push the latch open (fig. 44). This opening can be found on the closures on the femoral. panel and on the pelvic panels.

Comfort / Hygiene

Do not wear the orthosis directly on the skin, but instead over thin clothing (e.g. a t-shirt or a vest). Depending on the doctor's prescription it may be necessary to wear underwear over the orthosis, in case the orthosis may not be taken off during visits to the bathroom.

Material

Bandage:

48% Polyester 36% CoolMax® 16% Spandex

Orthosis component:

Steel splint PAG PF

Product Care / Cleaning / Maintenance

We recommend the use of a spare bandage (REF 23.601.XXX.XX) when the DvnaCox evolution has to be worn for an extended period of time.

Hand wash only.

Do not bleach.

Do not tumble dry.

Do not iron.

Do not dry-clean.

Storage / Disposal

Store the product in its original packaging in a dry place. Protect from dust, dirt and moisture. This product does not contain any environmentally hazardous substances. It can be disposed of with regular household waste.

CE-Conformity

This medical device complies with the requirements of the MDR, Regulation (EU) 2017/745.

Legend to Symbols

Manufacturer

Date of Manufacture

Distributor

Importer

EC Authorized Representative

Lot Number

Part Number

Serial Number

UDI Data Matrix

Follow User Manual

Caution

CE-Marking According to European Union Council Regulation 2017/745 Concerning Medical Devices (MDR)

Single Patient, for Multiple Use

Store in a Dry Place, Protect from Moisture

Protect from Sunlight

Temperature Limit

Humidity Limit

Atmospheric Pressure Limit, for Transport and Storage

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Wilhelm Julius Teufel GmbH Robert-Bosch-Straße 15 73117 Wangen Deutschland/Germany