Pavlik

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◆ Notizen

◆ Content

Description	EN-2
Indication	EN-2
Contraindication	EN-2
Function	EN-2
Reimbursement Code	EN-2
Medical Device	EN-2
Safety	EN-3
REF/Sizes	EN-4
Scope of Delivery	EN-4
Application/Adjustment	EN-5
Comfort/Hygiene	EN-7
Material	EN-7
Product Care/Cleaning/Maintenance	EN-7
Storage/Disposal	EN-8
Legend to Symbols	EN-8

◆ Description

With the Pavlik luxation bandage, the hip-knee bone flexion can be set to an angle of between 90 ° to 120 °. The bandage consists of a shoulder strap and a chest strap, and medial and lateral running leg straps with foot sections. Adjustments to the straps can be adapted according to the therapy course.

◆ Indication

- Hip-dysplasia
- Functional treatment of hip luxation

◆ Contraindication

- Non-compliant patients

◆ Function

- Dynamic functional method

◆ Reimbursement Code

DE HMV-No.: 05.05.01.1009

AT Tarif-No. 35012 012

◆ Medical Device

The Pavlik ist a medical device.

◆ ⚠ Safety

- Before this device is used, the orthosis has to be individually adapted to the patient by an experienced orthopedic technician and the patient has to be trained in the orthosis' use according to the doctor's instructions.
- Incorrect selection, adaptation, application, and/or failure to monitor the correct fit of this product may cause health damage.
- Check orthosis for damage (such as loose connections, fasteners, belts, etc.) each time it is put on.
- Do not remove the textile label from textile products.
- In case you suspect any performance malfunctions (e.g. an apparent reduction or lack of effect), immediately call on your doctor or orthopedic technician.
- All serious incidents that occur in connection with this product have to be reported to the Wilhelm Julius Teufel GmbH and the competent authority of the member state in which you are established. A „serious incident“ is any incident that caused, may have caused, or may cause, directly or indirectly, one of the following consequences:
 - a) the death of a patient, user, or other persons
 - b) temporary or permanent health deterioration of a patient, user or other persons
 - c) a serious public health threat.
- Never wear the orthosis on broken skin.
- The orthosis must be checked for proper function, correct fit and potential damage by a doctor or orthopedic technician at regular intervals. The timing of these intervals has to be defined individually for each patient.
- Should the product have been subject to disproportionate stress, it must be checked for potential damage by an orthopedic technician before further use.
- Harmful pressure caused by the orthosis may cause skin injuries or disrupt blood flow and must be avoided.
-  This product is intended for use on one person. It may not be reused on another person.

◆ REF/Sizes

REF	Size	Chest circumference	Foot size max.
23.400.011.00	S	30 – 40 cm	6 cm
23.400.012.00	S+	30 – 40 cm	7.5 cm
23.400.022.00	M	40 – 50 cm	7.5 cm

◆ Scope of Delivery

REF	Product
23.400.0XX.00	Pavlik Harness
23.402.000.00	Pavlik Harness Pad Kit
90.395.015.01	Pavlik_UserManual_DE_EN_FR_ES_IT

◆ Application/Adjustment

1. Fitting the Chest Section

1.1 Before initial use, the five buckles on the chest, shoulder and front leg belts should be opened. Lay the Pavlik Bandage open and stretched out on a suitable base (a nappy-changing table). The child should be positioned on its back, lying on the Pavlik Bandage so that the chest strap can be closed in a circular manner, a finger's width under the armpits whilst still above the nipples.

1.2 The shoulder straps should be crossed across the back so that the bandage does not slip over the shoulders.

2. Fitting the Leg Section

2.1 The leg straps are placed so that the longitudinal straps run parallel down both sides of the calves and over the inner and outer ankle. The belt part running underneath the foot must be positioned some centimetres in front of the rear edge of the heel, in the direction of the centre of the foot.

2.2 The leg straps running down the inner side are attached to the chest strap, so that a spreading of the upper thigh within the hip joints is possible to a degree of between 50° and 55° per leg.

Important: The attachment points of the leg straps must lie below an imaginary line drawn between the armpits. Following this, the outer leg straps are fastened to the rear buckles on the chest strap.

2.3 Should an increase in movement of the upper thighs be undesirable during the therapy course, the adjustable abduction inhibitor belt can be used. It has an insert sleeve for exchangeable stages, with the help of which an undesired adduction of the hips can be limited. A shorter stage is already enclosed within the belt, a longer stage is enclosed with the Pavlik bandage. The stages can also be shortened when required. To guarantee a correct fit of the foot parts, the three fasteners for the instep, above the heel and around the calf should not be applied too loosely.

2.3

3. Final fixing of the straps and padded straps

3.1 All belts can be shortened using a pair of scissors. However, make sure to leave enough excess end. Keep in mind that the belts might need to be changed during the course of the therapy. After the final length of the belts is set, the ends with the hook and loop fasteners of the padded flaps must also be fixed.

3.2 Threading the belt during fitting.

3.3 Threading the end of the belt.

3.3

3.4 The belt end is positioned as a deposit of the buckle and is shortened. With the hook and loop fastener the padded flap on the belt end can be finally fixed.

3.4

3.5 Correctly threaded fasteners with underlying padded straps. The tube-shaped pads of the shoulder belts can be pulled as far forward and downward as possible over the buckles.

3.5

When required, the pads on the shoulder and leg straps with the REF 23.402.000.00 can be reordered.

◆ Comfort/Hygiene

Do not wear the orthosis on broken skin. Wear the immobilizer over thin clothing.

◆ Material

Straps: PA velour

Pads: Velour foam

Foot sections: Subortholen, PA velour

◆ Product Care/ Cleaning/Maintenance

Hand wash only.

Do not bleach.

Do not tumble dry.

Do not iron.

Do not dry-clean.

◆ Storage/Disposal

This product does not contain any environmentally damaging materials. It can be disposed of as part of the normal household refuse.

◆ Legend to Symbols

Manufacturer

Date of Manufacture

Distributor

Importer

EC Authorized Representative

Lot Number

Part Number

Serial Number

UDI Data Matrix

Follow User Manual

Caution

CE-Marking According to European Union Council Regulation 2017/745 Concerning Medical Devices (MDR)

Single Patient,
for Multiple Use

Store in a Dry Place,
Protect from Moisture

Protect from Sunlight

Temperature Limit

Humidity Limit

Atmospheric Pressure
Limit, for Transport
and Storage

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