SACH-Plus Prothesenfuß

SACH Plus Prosthetic Foot

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User Manual

◆ Inhalt

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Contents

Description	EN-2
Technical Data	EN-3
Indication	EN-3
Medical Device	EN-3
Safety	EN-3
REF/Sizes	EN-5
Scope of Delivery	EN-5
Application/Adjustment	EN-5
Final Assembly	EN-6
Material	EN-6
Product Care/Cleaning/Maintenance	EN-6
Storage/Disposal	EN-6
Liability	EN-6
Warranty	EN-6
Legend to Symbols	EN-7

Description

The SACH Plus Foot offers exceptional attributes for a foot of its activity level. It is lighter than comparable feet since the titanium pyramid has been directly incorporated into the design. An additional pyramid adapter is not required. The SACH Plus Foot features an integrated toe spring made from composite material. The foot is available with three different toe resistances from which you can choose upon ordering. Heel strike can be adapted individually for each user, and three heel inserts of different firmness, which can be switched by the orthopedic technician, are included with each foot.

Technical Data

Weight of Foot/g	490 (Sz. 26, Pyramid Included)
Weight Limit/kg	115
Available Foot Sizes/cm	22 - 31
System Height/mm	88,3
Heel Height/mm	10
Toe Resistance Categories	2 - 4
Heel Insert Densities	Soft, Medium, Firm

Indication

Activity level: 1 + 2

Weight limit: 115 kg The user's adjusted body weight should not exceed the weight limit.

Approved for use in fresh water.

Adjusted body weight: The sum of the user's body weight and the weight of loads he normally or routinely carries. (See catalog for classification!)

The product is designed to be part of an individual prosthesis for a transfemoral. transtibial or knee-disarticulated amputee.

Medical Device

The SACH Plus Foot is a medical device.

- Before this component is used, both the component and the prosthesis it is used for have to be individually adapted to the user.
- Before the component is used, the user has to be trained in the component's use by an experienced orthopedic technician.
- Incorrect selection, adaptation, application and/or faulty inspection of this product may cause health damage.
- The prosthesis and its functional components must be checked for proper function, possible wear and tear and potential damage by an orthopedic technician at regular intervals.
- If the prosthesis or any of its components have been subject to disproportionate stress (e.g. due to a fall), they must be checked for potential damage by an orthopedic technician before further use.
- All serious incidents that occur in connection with this product have to be reported to the Wilhelm Julius Teufel GmbH and the competent authority of the member state in which you are established. A "serious incident" is any incident that caused, may have caused, or may cause, directly or indirectly, one of the following consequences: a) death of a patient, user, or other persons

- b) temporary or permanent health deterioration of a patient, user or other persons
- c) a serious public health threat.
- When using modular prosthetic components to fabricate a prosthetic limb, the following aspects have to be considered:
- Use components only according to their intended purpose.
- If components with different weight limits are combined, the weight limit of the weakest component applies to the whole prosthesis.
- If components for different activity levels are combined, the activity level of the component with the lowest activity level applies to the whole prosthesis.
- Using individual, certified components with the CE-marking does not release technicians from their duty to check component combinations within their means for fitness for purpose, proper assembly and safety.
- Should there be evidence that a certain combination of components is not as safe as required, the components must not be combined.
- The prosthesis has to be configured and assembled according to the recognized technical rules of the trade.
- The locking torque specifications given for the modular component have to be met. A suitable tool (i.e. a torque wrench) has to be used to set the locking torque accordingly.
- Before the prothesis is handed over to the user, all modular components have to be secured against accidental

- loosening or twisting with suitable threadlocking fluid.
- Safety related regulations for individual components (e.g. specific combination possibilities, maintenance intervals etc.) must be followed.
- If different components (from one or from several manufacturers) are combined: Only use components which meet the requirements of DIN EN ISO 10328 and the Medical Devices Act.
- In the event of damage to a prosthesis which combines components from different manufacturers, every manufacturer can only be held liable for the failure of those components produced by them. Beyond that, manufacturers can only be held liable if their components are proven to have caused damage to or failure of components which were produced by other manufacturers.
- This product is intended for use on one person. It may not be reused on on another person.

RFF/Sizes

Please see the current edition of the Teufel "Prosthetics Catalog".

Scope of Delivery

REF	Product
WW.SFP-XXXX-XX	SACH Plus Foot
90.395.101.00	SACHPlus_UserMa- nual_DE_EN

Application/ Adjustment

We recommend following the assembly instructions in the order suggested below.

Assembly

Follow the fundamental assembly rules

- · Consider the intended heel height of the prosthesis
- Consider the position of the limb and the configuration of the prosthesis that is based on it

Static Assembly

Recommended alignment reference line: TKA (Trochanter- Knee - Ankle)

Installing and Switching the Heel Inserts

The SACH Foot comes with three heel inserts of different firmness. It is thus possible to adjust the heel firmness individually.

Tip: Start by installing the medium firm heel insert. This will usually result in having to switch the insert no more than once. To install, place the insert in the cut-out on the heel, it will automatically lock in place. To remove, use a flat screwdriver to slightly lift the heel insert out of its recess. The insert can then be easily pulled out and removed by hand.

- if the medium heel insert is too soft while test-walking; swap medium insert for hard insert
- if the medium heel insert is too firm while test-walking; swap medium insert for softer insert

Incl. 3 Heel Inserts

Final Assembly

The screw connection between the foot pyramid and the selected tube clamp adapter or pylon adapter has to be fastened with the torque required for said adapter.

Secure all screws with medium strength threadlocking fluid for the definite final assembly.

Material

Titanium Composite Stainless steel Plastic

Product Care/ Cleaning/Maintenance

If necessary, SACH Feet can be cleaned with a damp cloth.

Conduct an annual visual inspection and a performance test of the foot.

If required, send the foot to us for maintenance.

Storage/Disposal

Store the product in its original packaging in a dry place. Protect from dust, dirt and moisture. This product does not contain any environmentally hazardous

substances. It can be disposed of with household refuse.

Liability

The manufacturer recommends using the product only under the specified conditions and for the intended purposes, as well as only using modular components which have been certified for the prosthesis, and to maintain the product according to the user manual. The manufacturer will not accept liability for damages caused through improper use of the components.

Warranty

The warranty period for SACH Plus Feet is 12 months.

Legend to Symbols

Manufacturer

Date of Manufacture

Distributor

Importer

EC Authorized Representative

Lot Number

Part Number

Serial Number

UDI Data Matrix

Follow User Manual

Caution

CE-Marking According to European Union Council Regulation 2017/745 Concerning Medical Devices (MDR)

Single Patient, for Multiple Use

Store in a Dry Place, Protect from Moisture

Protect from Sunlight

Temperature Limit

Humidity Limit

Atmospheric Pressure Limit, for Transport and Storage

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90.395.101.00 Rev. 2022/08 © Wilhelm Julius Teufel GmbH
Wilhelm Julius Teufel GmbH Robert-Bosch-Straße 15 73117 Wangen Deutschland/Germany Phone: +49 (0)7161 15684-0
WillowWood Global LLC 15441 Scioto Darby Road Mt. Sterling, OH 43143, US