Modularadapter Modular Components

Gebrauchsanweisung User Manual

Produktbeispiel | Product Example

www.teufel-international.com

◆ Inhalt

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◆ Notizen

◆ Contents

Description	8
Specification	8
Indication	9
Contraindication	9
Safety	10
Legend to Symbols	13

◆ Description

Modular prosthetic components for the fabrication of prosthetic limbs.

◆ Specification

Weight Limit:

The weight limit is indicated on the product label and in our prosthetic catalogue.

Activity Levels:

Titanium adapters are approved for activity levels 1-4.

Stainless steel adapters are approved for activity levels 1-4.

Aluminum adapters are approved for activity levels 1-3.

Please also note the explanation of activity levels and weight limits in our prosthetics catalogue, as well as the technical data (i.e. locking torques) the adapters listed therein.

◆ Indication

- Modular prosthetic limbs.
- Waterproof prosthetic limbs, if all used components are marked waterproof.

◆ Contraindication

- Use of components for non-approved activity levels or weight limits.
- Use of non-waterproof components in waterproof prostheses.

Safety

- When using modular prosthetic components to fabricate a prosthetic limb, the following aspects have to be considered:
- Use components only according to their intended purpose.
- If components with different weight limits are combined, the weight limit of the weakest component applies to the whole prosthesis.
- If components for different activity levels are combined, the activity level of the component with the lowest activity level applies to the whole prosthesis.
- Using individual, certified components with the CE-marking does not release technicians from their duty to check component combinations within their means for fitness of purpose, proper assembly and safety.
- Should there be evidence that a certain combination of components is not as safe as required, the components must not be combined.
- The prosthesis has to be configured according to the recognized technical rules of the trade.
- The locking torque specifications given for the modular component have to be met. A suitable tool (i.e. a torque wrench) has to be used to set the locking torque accordingly.
- The relevant locking torques have either been

lasered directly onto the adapter, or can be found in the prosthetic catalogue.

- If the locking torque has been lasered directly onto the adapter, the following aspects have to be considered:
- If different locking torques have to be considered for different screws on the adapter, the relevant torque is indicated on each individual screw.
- If the same locking torque applies to different screws on the adapter, the relevant torque is only given once on the adapter.
- Locking torques always apply to the screw threads of the adapters. If two adapters are connected with each other, the rele-

vant locking torque is the one given for the adapter that has a screw thread. Because of this, some adapters may not have any applicable locking torque at all, since they feature no screw threads (e.g. 4-hole-adapter with pyramid).

- Before the prosthesis is handed over to the user, all modular components have to be secured against accidental loosening or twisting with suitable thread-locking fluid.
- Safety related regulations for individual components (e.g. specific combination possibilities, maintenance intervals etc.) must be followed.

- If different components (from one or from several manufacturers) are combined: Only use components which meet the requirements of DIN EN ISO 10328 and the Medical Devices Act.
- In the event of damage to a prosthesis which combines components from different manufacturers, every manufacturer can only be held liable for the failure of those components produced by them.

Beyond that, manufacturers can only be held liable if their components can be proven to have caused damage to or failure of components which were produced by other manufacturers.

- The modular components may not be reused on other persons.

◆ Legend to Symbols

Manufacturer

EC Authorized
Representative

Lot Code

Part Number

Serial Number

Follow User Manual

Caution

CE-Marking According
to Council Directive
93/42/EEC Concerning
Medical Devices

◆ Notes

 **Notizen/Notes**

<hr/>	<hr/>
<hr/>	<hr/>
<hr/>	<hr/>
<hr/>	<hr/>
<hr/>	<hr/>
<hr/>	<hr/>
<hr/>	<hr/>
<hr/>	<hr/>
<hr/>	<hr/>
<hr/>	<hr/>
<hr/>	<hr/>
<hr/>	<hr/>
<hr/>	<hr/>
<hr/>	<hr/>
<hr/>	<hr/>
<hr/>	<hr/>
<hr/>	<hr/>
<hr/>	<hr/>
<hr/>	<hr/>
<hr/>	<hr/>
<hr/>	<hr/>
<hr/>	<hr/>
<hr/>	<hr/>
<hr/>	<hr/>

◆ Notizen/Notes

_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____

Wilhelm Julius Teufel GmbH
Robert-Bosch-Straße 15
73117 Wangen
Deutschland/Germany
Phone: +49 (0)7161 15684-0
Fax: +49 (0)7161 15684-333
www.teufel-international.com