
# Modularadapter Modular Components

# Gebrauchsanweisung User Manual

Produktbeispiel | Product Example

[www.wjt-ortho.com](http://www.wjt-ortho.com)

## ◆ Inhalt

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## ◆ Contents

Description	EN-2
Technical Data	EN-2
Indication	EN-3
Contraindication	EN-3
Medical Device	EN-4
Safety	EN-4
Specific	
Safety Instructions	EN-7
Application/Adjustment	EN-8
Alignment	EN-8
Final Assembly	EN-10
Materials	EN-11
Product Care/ Cleaning/Maintenance	EN-11
Storage/Disposal	EN-11
Liability	EN-12
Warranty	EN-12
Legend to Symbols	EN-12
ISO-Standard	EN-13

## ◆ Description

Modular prosthetic components for the fabrication of prosthetic limbs.

## ◆ Technical Data

### **Weight Limit:**

The weight limit is indicated on the product label and in our prosthetic catalog.

**Activity Levels:**

Titanium adapters are approved for activity levels 1-4.

Stainless steel adapters are approved for activity levels 1-4.

Aluminum adapters are approved for activity levels 1-3.

Please also note the explanation of activity levels and weight limits in our prosthetics catalog, as well as the technical data (i.e. locking torques) for the adapters listed therein.

**◆ Indication**

- Modular prosthetic limbs.
- Waterproof prosthetic limbs, if all used components are marked as waterproof.

**◆ Contraindication**

- Use of components for non-approved activity levels or weight limits.
- Use of non-waterproof components in waterproof prostheses.

## ◆ Medical Device

The modular components are a medical device.

## ◆ ⚠ Safety

- Before this component is used, both the component and the prosthesis it is used for have to be individually adapted to the user.
- Before the component is used, the user has to be trained in the component's use by an experienced orthopedic technician.
- Incorrect selection, adaptation, application and/or faulty inspection of this

product may cause health damage.

- The prosthesis and its functional components must be checked for proper function, possible wear and tear and potential damage by an orthopedic technician at regular intervals.
- If the prosthesis or any of its components have been subject to disproportionate stress (e.g. due to a fall), they must be checked for potential damage by an orthopedic technician before further use.
- All serious incidents that occur in connection with this product have to be

reported to the Wilhelm Julius Teufel GmbH and the competent authority of the member state in which you are established. A "serious incident" is any incident that caused, may have caused, or may cause, directly or indirectly, one of the following consequences:

- a) death of a patient, user, or other persons
  - b) temporary or permanent health deterioration of a patient, user or other persons
  - c) a serious public health threat.
- When using modular prosthetic components to fabricate a prosthet-

ic limb, the following aspects have to be considered:

- Use components only according to their intended purpose.
- If components with different weight limits are combined, the weight limit of the weakest component applies to the whole prosthesis.
- If components for different activity levels are combined, the activity level of the component with the lowest activity level applies to the whole prosthesis.
- Using individual, certified components with the CE-marking does not

release technicians from their duty to check component combinations within their means for fitness for purpose, proper assembly and safety.

- Should there be evidence that a certain combination of components is not as safe as required, the components must not be combined.
- The prosthesis has to be configured and assembled according to the recognized technical rules of the trade.
- The locking torque specifications given for the modular component have to be met. A suitable tool (i.e. a torque wrench) has to be

used to set the locking torque accordingly.

- Before the prosthesis is handed over to the user, all modular components have to be secured against accidental loosening or twisting with suitable Threadlocking fluid.
- Safety related regulations for individual components (e.g. specific combination possibilities, maintenance intervals etc.) must be followed.
- If different components (from one or from several manufacturers) are combined: Only use components which meet the require-

ments of DIN EN ISO 10328 and the Medical Devices Act.

- In the event of damage to a prosthesis which combines components from different manufacturers, every manufacturer can only be held liable for the failure of those components produced by them. Beyond that, manufacturers can only be held liable if their components are proven to have caused damage to or failure of components which were produced by other manufacturers.
-  This product is intended for use on one person. It may not be reused on another person.

## Specific Safety Instructions

- The relevant locking torques have either been lasered directly onto the adapter, or can be found in the prosthetic catalogue.
- If the locking torque has been lasered directly onto the adapter, the following aspects have to be considered:
- If different locking torques have to be considered for different screws on the adapter, the relevant torque is indicated on each individual screw.
- If the same locking torque applies to different screws on the adapter,

the relevant torque is only given once on the adapter.

- Locking torques always apply to the screw threads of the adapters. If two adapters are connected with each other, the relevant locking torque is the one given for the adapter that has a screw thread. Because of this, some adapters may not have any applicable locking torque at all, since they feature no screw threads (e.g. 4-hole-adapter with pyramid).

## ◆ Application/ Adjustment

⚠ Follow the instructions for application and adjustment below.

## ◆ Alignment

⚠ Follow the recognized rules of the trade during the alignment of the prosthesis.

⚠ Follow the specifications on the weight limit of the individual adapters. Keep in mind that the component with the lowest weight limit defines the weight limit for the whole prosthesis.

 Follow the instructions regarding maximum locking

torque as indicated on the adapters.

REF	Screws	Torque/Nm
<b>Ww.700-200</b>	M6x16	12
<b>Ww.700-250</b>	M6x16	12
<b>Ww.700-300</b>	M6x20	12
<b>Ww.700-375</b>	M6x25	12
<b>Ww.FND-135105</b>	M6x16 M6x10	12 12
<b>Ww.FND-135107</b>	M6x16 M6x10	12 12
<b>Ww.FND-136124</b>	M8x16	15
<b>Ww.FND-136146</b>	M8x16	15
<b>Ww.FND-224120</b>	M6x16	12
<b>Ww.FND-224150</b>	M6x16 (Akt. 1-3) M6x35 (Akt. 4)	12 12

REF	Screws	Torque/Nm
WW.MHD-DE-P	-	-
WW.PPF-1208	M8x18	20
WW.PPF-1209	M6x12	4
WW.PPF-FS006	M8x18	20
WW.PPF-GKIT	M8x55	20
WW.PPF-700250	M5x12	7
WW.PPF-994114	M5x20	7

## ◆ Final Assembly

⚠ All screwed connections have to be secured with the indicated locking torques. Use medium-strength threadlocking fluid to ensure maximum security.

⚠ When installing a pylon adapter on the tube clamp adapter, make sure to de-grease the clamping surfaces of both adapters

before their final connection. When securing the connection with the correct torque, this will ensure that the pylon adapter will not accidentally twist inside the tube clamp adapter.

## ◆ Materials

Titanium  
Stainless steel  
Aluminum

## ◆ Product Care/ Cleaning/ Maintenance

If necessary, modular adapters can be cleaned with a damp cloth.

Conduct an annual visual inspection and performance check of the prosthesis and the modular adapters that were used to build it.

## ◆ Storage/ Disposal

Store the product in its original packaging in a dry place. Protect from dust, dirt and moisture. This product does not contain any environmentally hazardous substances. It can be disposed of with regular household waste.

## ◆ Liability

In the event of damage, every manufacturer can only be held liable for the components produced by them. Any extended liability is only possible if the components from the manufacturer in question can be proven to have caused damage to or loss of function of the components produced by another manufacturer.

## ◆ Warranty

The warranty for modular components is 12 months.

## ◆ Legend to Symbols

Manufacturer

Date of Manufacture

Distributor

Importer

EC Authorized Representative

Lot Number

Part Number

Serial Number

UDI Data Matrix

Follow User Manual

Caution

CE-Marking According  
to European Union  
Council Regulation  
2017/745 Concerning  
Medical Devices (MDR)

Single Patient,  
for Multiple Use

Store in a Dry Place,  
Protect from Moisture

Protect from  
Sunlight

Temperature Limit

Humidity Limit

Atmospheric  
Pressure Limit, for  
Transport and Storage

Medical Device

## ◆ ISO-Standard

The modular adapters have  
been tested by a certified  
testing institute according  
to DIN EN ISO 10328.

Importeur:  
Importer:

Wilhelm Julius Teufel GmbH  
Robert-Bosch-Straße 15  
73117 Wangen  
Deutschland/Germany  
Phone: +49 (0)7161 15684-0  
Fax: +49 (0)7161 15684-333  
[www.wjt-ortho.com](http://www.wjt-ortho.com)

#### WILLOWWOOD.

WillowWood Global LLC  
15441 Scioto Darby Road  
P.O. Box 130  
Mt. Sterling, Ohio 43143  
USA

[www.willowwood.com](http://www.willowwood.com)

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© Wilhelm Julius Teufel GmbH