Sensor Knee

Gebrauchsanweisung User Manual

◆ Inhalt

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
◆ Technical Data

1. Frame
2. Hydraulic Unit
3. Knee Ball with Pyramid Receiver
4. Tube Clamp Attachment, with s4 Hexagon Socket Screw, 10 Nm
5. Silver Flexion Resistance Adjustment Knob
6. M6 Hexagon Socket Screw
7. Black Yield Resistance Adjustment Knob
8. Yield/Lock Lever
9. s2 Hexagon Socket Screw, for Adjustment of End Position Damping
10. Plastic Insert
11. Axle Bolt (Alignment Reference Point for Knee Pivot Point)
12. s4 Set Screws, 10 Nm
13. Oil Level Gauge (Golden Ball End)

REF	Product	Total Length /mm	Weight /g	Maximum Flexion	Weight Limit /kg	Construction Height of Knee Pivot Point with Reference to Socket Attachment Block
OE.2SR500	Sensor Knee Joint Silver	276	1141	120°	125	40 mm
OE.2SR502-BLK	Sensor Knee Joint Black	276	1141	120°	125	40 mm

◆ Indication

Activity Level: 2 – 3+

Weight Limit: 125 kg

The weight limit is the sum of the user's body weight, the weight of the user's clothes, and the weight of loads he normally or routinely carries. (See Prosthetics Catalog for classifications.)

The product is designed to be part of an individual prosthesis for one transfemoral or knee-disarticulated amputee.

◆ Medical Device

The Sensor Knee Joint is a medical device.

◆ ⚠ Safety

- Before this component is used, both the component and the prosthesis it is used for have to be individually adapted to the user.
- Before the component is used, the user has to be trained in the component's use by an experienced orthopedic technician.
- Incorrect selection, adaptation, application and/or faulty inspection of this product may cause health damage.
- The prosthesis and its functional components must be checked for proper function, possible wear and tear and potential damage by an orthopedic technician at regular intervals.
- If the prosthesis or any of its components have been subject to disproportionate stress (e.g. due to a fall), they must be checked for potential damage by an orthopedic technician before further use.
- All serious incidents that occur in connection with this product have to be reported to the Wilhelm Julius Teufel GmbH and the competent authority of the member state in which you are established. A "serious incident" is any incident that caused, may have caused, or may cause, directly or indirectly, one of the following consequences:
 - a) death of a patient, user, or other persons
 - b) temporary or permanent health deterioration of a patient, user or other persons
 - c) a serious public health threat.
- When using modular prosthetic components to fabricate a prosthetic

limb, the following aspects have to be considered:

- Use components only according to their intended purpose.
 - If components with different weight limits are combined, the weight limit of the weakest component applies to the whole prosthesis.
 - If components for different activity levels are combined, the activity level of the component with the lowest activity level applies to the whole prosthesis.
 - Using individual, certified components with the CE-marking does not release technicians from their duty to check component combinations within their means for fitness for purpose, proper assembly and safety.
 - Should there be evidence that a certain combination of components is not as safe as required, the components must not be combined.
 - The prosthesis has to be configured and assembled according to the recognized technical rules of the trade.
 - The locking torque specifications given for the modular component have to be met. A suitable tool (i.e. a torque wrench) has to be used to set the locking torque accordingly.
 - Before the prosthesis is handed over to the user, all modular components have to be secured against accidental loosening or twisting with suitable Threadlocking fluid.
 - Safety related regulations for individual components (e.g. specific combination possibilities, maintenance intervals etc.) must be followed.
- If different components (from one or from several manufacturers) are combined: Only use components which meet the requirements of DIN EN ISO 10328 and the Medical Devices Act.
 - In the event of damage to a prosthesis which combines components from different manufacturers, every manufacturer can only be held liable for the failure of those components produced by them. Beyond that, manufacturers can only be held liable if their components are proven to have caused damage to or failure of components which were produced by other manufacturers.
 -  This product is intended for use on one person. It may not be reused on another person.

◆ REF/Sizes

REF	Product
OE.2SR500	Sensor Knee Joint Silver
OE.2SR502-BLK	Sensor Knee Joint Black

◆ Scope of Delivery

REF	Product
OE.2SR500 / OE.2SR502-BLK	Sensor Knee Joint Silver Sensor Knee Joint Black
90.395.128.00	User Manual

◆ Accessories


• Optional Accessories

REF	Product
21.989.003.00	Foam Cover AK M3
AP.TI-110LONG	Pylon Connector Pyr.Receiver Height 19 166kg Ti
AP.TI-112-P-W	Pylon Adap. w/ Pyr. 88g166kgØ30mm Titan
AP.TI-115DM-W	Dual Adapter Pyramid Height 22 166kg Ti Wp
WW.MHD-DE-P	Dual Pyramid Adapter Height 15 166kg Ti

◆ Application/Adjustment

• Recommended Prosthesis Alignment

1. The alignment reference line should run directly through the knee's pivot point (fig. 2).
2. The knee's pivot point is positioned 8 mm behind and 5 mm above the joint axis (axle bolt) (alignment reference line for the knee's pivot point, fig. 1, section 11).

3. The load line should also run through the knee's pivot point.
4. If the user is less active, we recommend changing the load line so it runs through a point a little anterior to the knee's pivot point. If the user is more active, change the load line so it runs through a point a little posterior to the knee's pivot point.

• Static Fitting

Upon delivery, every joint is set to a basic setting. Ask the user to stand up. Remain close, so you can support the user if necessary. Individually adjust all components until the user is able to stand securely. Possible

adjustments are described in the section titled "Dynamic Fitting".

To start the dynamic fitting process, first remove the black cap by pulling it off towards posterior (fig. 3).

• Dynamic Fitting

Yield/Lock Lever:

Stance Phase Control

Users can choose from 2 stance phase modes:

Lock Mode: The lock mode is activated by the user's body weight, it secures the joint for a flexion angle of up to 15-20°. Once the joint is unloaded, swing phase is initiated and the lock gets deactivated.

Yield Mode: During stance phase, the joint yields with controlled, user-specific resistance.

 **Practical Tip:**

Users can easily switch between the two modes. When users first start wearing a prosthesis with a Sensor Knee, we recommend the following basic settings:

Lock Mode: For users who prefer a high level of stability and security.

Yield Mode: For more active and secure walkers, or users who are active enough to climb and descend stairs step-over-step and to walk on slopes.

Activating Lock Mode (Fig. 4)

Push the yield/lock lever **down**
= lock setting (fig. 4)

- ideal for uneven terrains, provides security
- recommended for older or hip-disarticulated amputees

Activating Yield Mode (Fig. 5)

Push the yield/lock lever **up**
= Yield mode

- Ideal for active users and secure walkers.
- Makes it possible to climb and descend stairs step-over-step, and allows users to walk on slopes.

Adjusting the Hydraulic Unit

The hydraulic unit of the Sensor Knee offers four optional settings:

- Response threshold (s4 hexagonal socket screw, fig. 1, no. 6)
- Yield resistance (black adjustment knob, fig. 1, no. 7)
- Flexion resistance (silver adjustment knob, fig. 1, no. 5)
- End position damping (s2 hexagon socket screw at the distal end of the tube clamp, fig. 1, no. 9)

Adjusting Sensitivity (Fig. 6)

- The response behavior of the stance phase control (sensitivity) can be adjusted with the s4 hexagon socket screw. This function is important for both stance phase modes.
- For easier handling we recommend using a hexagonal wrench with a dome head.
- The sensitivity setting determines how much weight it takes to activate stance phase control upon heel strike.
-  **Technician's Tip:** In case the joint "jams" at toe-off, the sensitivity of the joint has to be decreased.
- Perform a heel strike test with the user while the joint is in lock mode in order to check whether the sensitivity setting is adequate.
- **Turning the Screw Clockwise -** Increases weight sensitivity - for users with a higher body weight or users with high ground reaction forces.
- **Turning the Screw Counterclockwise -** Decreases weight sensitivity - for users with a lower body weight or users with low ground reaction forces.


Adjusting Yield Resistance - Black Adjustment Knob (Fig.7)

- Controls the dynamic stance phase damping for knee flexion while the user walks on stairs or on slopes.

Flip the yield/lock lever to the "up"- position. Adjust the dynamic stance phase damping with the black adjustment knob (fig.7):

- **Turning the Knob Clockwise -** Increases resistance.
- **Turning the Knob Counterclockwise -** Decreases resistance.

 **Caution!** The adjustment knob must never be tightened all the way, since this might damage the hydraulic unit. To check this, carefully completely tighten the knob, then loosen it again by at least 10°.


Adjusting Flexion Resistance - Silver Adjustment Knob (Fig. 8)

- Limits flexion, i.e. the heel amplitude, of the lower leg during swing phase.

Adjust the flexion resistance as desired with the silver adjustment knob (fig.8):

- **Turning the Knob Clockwise -**
Decreases the maximum flexion angle.
- **Turning the Knob Counterclockwise -**
Increases the maximum flexion angle.

⚠ Caution! The adjustment knob must never be tightened all the way, since this might damage the hydraulic unit. To check this, carefully completely tighten the knob, then loosen it again by at least 10°.


Adjusting the End Position Damping - s2 Hexagon Socket Screw (Fig. 9)

- The sensor knee progressively dampens the last 5° of extension. This helps prevent an abrupt extension stop without limiting motion during swing phase.
- The basic factory settings have to be changed only very rarely.

To change the settings of the end position damping, pull the pylon adapter off the joint. You will now see the screw with which you can adjust the end position damping.

- **Turning the Screw Clockwise -**
Increases the end position damping
- **Turning the Screw Counterclockwise -**
Decreases the end position damping
- The user benefits from the lower effort it takes to extend the knee.

When the user is walking on stairs and on slopes the following has to be considered:

- Yield mode is recommended (yield/lock lever at upper position)
- Allow the user to practice putting weight on the joint and to develop a feeling for the yield function.
- Adjust resistance, if necessary.
- Have the user start with the bottom step when he first practices walking on stairs.
- During the final adjustments of the joint settings the user should not be using any crutches.
- Make sure the user is sufficiently secured during the final adjustment.
- Weight has to be put on the joint to activate the dynamic stance phase damping.
- While practicing, heel strike should happen while the joint is extended, to ensure that force is transmitted properly into the joint.
- You have adjusted the damping resistance correctly when the yield speed of the healthy leg and the prosthetic knee are the same when the user is descending stairs.

When the user is sitting down while wearing the prosthesis, the following has to be considered:

The user can also benefit from the support of the Sensor Knee's hydraulic unit when sitting down.

- The user should place the prosthetic foot about 1/4 of the length of a foot in front of the healthy foot.
- When putting weight on the prosthesis, the yield function will now support the user when sitting down.

◆ Final Assembly

Do not use talcum powder or any other powder. Powders bind fat and thus will cause the joint to run dry. When putting a cosmetic foam cover on the prosthesis, silicone spray may be used.

The following torque values should be applied:

Tube clamp screw, s4: **10 Nm**.

Set screws, s4: **10 Nm**.

Secure the screws with medium strength threadlocking adhesive (such as Loctite 243).

⚠ Caution! The set screws of the proximal pyramid receiver must never disappear in the frame, since otherwise they will not be held securely by the screw threads. In case the short set screws do disappear inside the frame, use the long set screws which are also included in delivery.

◆ Material

Titanium
Aluminum
Stainless Steel
Plastic

◆ Product Care/ Cleaning/Maintenance

Oil Level Check of the Hydraulic System

The golden ball head between the black and the silver adjustment knob indicates the oil level of the hydraulic system. Temperature and pressure changes will influence the position of the ball head. It may protrude or even be set back a little.

⚠ Caution! The ball head must never be significantly set back into the housing. If this is the case, the hydraulic system may be losing oil. The joint must not be used any more and should be sent to us for inspection.

If necessary, Sensor Knee joints can be cleaned with a damp cloth. Conduct an annual visual inspection as well as a performance test. We recommend sending the joint to us once a year for maintenance.

◆ Storage/Disposal

Store the product in its original packaging in a dry place. Protect from dust, dirt and moisture. This product does not contain any environmentally hazardous substances. It can be disposed of with regular household waste.

◆ Legend to Symbols

Manufacturer

Date of Manufacture

Distributor

Importer

EC Authorized Representative

Lot Number

Part Number

Serial Number

UDI Data Matrix

Follow User Manual

Caution

CE-Marking According to European Union Council Regulation 2017/745 Concerning Medical Devices (MDR)

Single Patient, for Multiple Use

Store in a Dry Place, Protect from Moisture

Protect from Sunlight

Temperature Limit

Humidity Limit

Atmospheric Pressure Limit, for Transport and Storage

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www.wjt-ortho.com

Ortho Europe Ltd,
Ability House, Nuffield Way
Abingdon, Oxfordshire, OX14 1RL, UK
www.ortho-europe.com

Assistive Technology Group
Unit 12, Deansgrange Business Park
Blackrock, Dublin, A94 HF60
Ireland

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