# ESP 3 Gänge Ventile

ESP Screw Valves

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User Manual

## ◆ Inhalt

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## ◆ Contents

Description	6
Quick Reference Guide	6
Accessories	8
Product Care/ Cleaning/Maintenance	8
Liability	8
CE-Conformity	8

## ◆ Description

ESP Screw Valves are designed to be used in AK prostheses.

## ◆ Quick Reference Guide

Scope of Delivery: Valve body (REF ES.P100), with required accessories, valve for selected version (REF ES.P1, ES.P2, ES.P3 or connector ES.P4 for elevated vacuum systems) including dummy and screw for attaching the dummy.

1. Positioning: Position the valve in the distal area of the socket (in a spot where the surface is not too curved).
2. Put a layer of stockinette (where applicable) and/or a PVA sleeve on the model.

3. Use a permanent marker to mark the center of the dummy hole. Reinforce the PVA sleeve with a piece of PE tape.
4. Position the dummy on the model in the desired spot.
5. Vacuum formed sockets: Position the dummy on the model in the desired spot. Screw it into place. Ensure sufficient air exhaustion, especially in the area around the dummy. If necessary, add additional drill holes to the model to ensure the efficient exhaustion of air.
6. Laminated sockets: Position the dummy on the model in the desired spot. Screw it into place.
7. Once the material has cooled/cured, sand all excess material off until the surface of the screw is exposed. Remove the screw. Continue to carefully sand the material down, working parallel to the surface of the dummy, until you have exposed the surface of the dummy.
8. Insert the valve body from the inside of the socket (including the black seal).
9. Take the valve body and the threaded ring and screw them together (left-hand thread!).
10. Now screw the valve into the valve body (right-hand thread!).
11. ESP Screw Valves are not waterproof.

## ◆ Accessories

The valves REF ES.P1, ES.P2, ES.P3 and the elevated vacuum system connector ES.P4 are available separately.

The valve body REF ES.100, which fits all versions of the valve, is also sold separately.

## ◆ Product Care/ Cleaning/ Maintenance

The performance of the ESP Screw Valves must be checked at regular intervals (depending on the intensity with which the prosthesis is used). Dirt inside the air channels may impair the functionality of the ESP Screw Valves. Remove the expulsion valve and clean it. Then re-insert the expulsion valve. If the screw valves should be damaged or worn-out, we recommend replacing them.

## ◆ Liability

The manufacturer recommends that the products be only used for the intended purpose. Also, only combinations of modular components which were tested for the prosthesis should be used.

The manufacturer will not accept liability for any damages caused by components which were not approved for this product by the manufacturer.

## ◆ CE-Conformity

This product fulfills the requirements of the Council Directive 93/42/EEC concerning medical devices.

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Wilhelm Julius Teufel GmbH
ESP, LLC.